

# Precision® Neonatal Foam Band Instructions



Center label on info area to allow best possible seal. Press firmly on label to secure.



Remove white backing from adhesive on shield.



Press shield down firmly from left to right on flat surface several times.



Apply foam band to patient.

**WARNING:** Foam band may restrict circulation if applied too tightly, resulting in nerve damage or tissue damage.

When banding the infant, the foam band should fit snugly (similar to a watch band). In this way, when the infant loses excess body fluid during the first 24 hours after birth\*, the band will not become too loose.

\*Reference: Maternity and Gynecologic Care. The Nurse and the Family Jensen Bobak 1985 (pg. 679)