

Spee-D-Cool™ Gel Pads



Reduce patient discomfort before and after a biopsy.

REUSABLE!

Original Spee-D-Cool™

SCP-30
3" Circle
48/Box

SCP-60
4" Quadrant
24/Box

Tropical Design


Ensure your patients have a positive experience while marketing your facility!

SCP-60
4" Quadrant
24/Box

Soft and pliable even when frozen

- Its soft rounded edge conforms naturally to a woman's anatomy
- Provides 30 minutes of soothing cold therapy
- The soft brushed fabric surface acts as a thermal buffer to eliminate the shock of initial contact — no need to protect the skin with a thick gauze pad
- Slim design allows pad to be worn under clothing without being noticed
- 3" circle includes latex-free adhesive strips for secure positioning
- Available in 2 sizes; 4" quadrant also available in tropical design and with custom imprint
- Latex-free

Gel Pads with Removable Sleeves



Our heart-shaped gel pads with sleeves offer a fresh cover as needed during patient recovery. Removable sleeves prolong the use of gel pads. Spee-D-Cool™ gel pads with sleeves offer the same functionality as the original gel pads, with the added benefit of extended reusability.

Gel pads can be placed within a soft brushed fabric sleeve, eliminating the shock of initial contact. Plastic gel pads can be easily cleaned with soap and water. Each gel pad comes with 2 fabric sleeves, one of each design.



SCPS-50
24 Heart-shaped gel pads/Box
48 sleeves/Box



SCPSLV1
Heart-shaped sleeve
ribbon design #1
24 sleeves/Pack



SCPSLV2
Heart-shaped sleeve
ribbon design #2
24 sleeves/Pack

Get Free Samples! Contact Us Today 800.435.4242
or visit pdchealthcare.com/mammo to view our complete line of Mammography products