

StainerShield®

Slide Labeling Solutions




Tired of spending valuable time and money on misidentified specimens, redraw samples and repeat testing? StainerShield® slide labels are a simplified solution to the complex problem of properly identifying and tracking specimen slides in a laboratory environment. Lab technicians can now print on-demand, patient-specific labels and apply them to tissue and specimen slides as the slides are produced, before they go through the stainer. No over-laminate flap, UV ink, or etching system is needed.


These labels can be printed on-demand using standard desktop thermal printers. In one simple printing step, technicians can print slide labels with both text and bar codes — including 2D bar codes — that contain data directly from the lab's information system.




**StainerShield® Slide Labels are
Xylene and Stain Resistant**



Product Number	Label Size	Carrier Width	Core Size	Wound	Max OD	Material/Adhesive	Qty.
TDSS1-1-7878	7/8" x 7/8"	1"	1"	In	4-13/32"	StainerShield® Direct Thermal/ Glove-Friendly Permanent	1,800/Roll- 4Rolls/Box
TDSS15-17878			1-1/2"		4-13/32"		
TDSS3-1-7878			3"		5"		



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TDSS1-2-7878	7/8" x 7/8"	2"	1"	In	StainerShield® Direct Thermal/Glove-Friendly Permanent	3,600/Roll- 2Rolls/Box
TDSS3-2-7878			3"			



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TDSS1-4-7878	7/8" x 7/8"	4"	1"	In	StainerShield® Direct Thermal/ Glove-Friendly Permanent	7,200/Roll- 1Roll/Box
TDSS3-4-7878	7/8" x 7/8"		3"			
TDSS1-4-1278	7/8" x 1/2"		1"			
TDSS3-4-1278	7/8" x 1/2"		3"			



Product Number	Label Size	Carrier Width	Core Size	Wound	Material/Adhesive	Qty.
TDSS1-4-1516	15/16" x 15/16"	3-7/8"	1"	In	StainerShield® Direct Thermal/Glove-Friendly Permanent	7,200/Roll- 1Roll/Box
TDSS3-4-1516			3"			

Contact Us Today 800.435.4242

We offer many more product categories and customization options.

EDI ordering via GHX (Global Healthcare Exchange)

PDC Healthcare 27770 N. Entertainment Dr. Suite 200 Valencia, CA 91355

Tel 800 435 4242 | Fax 800 321 4409 | info@pdchealthcare.com | pdchealthcare.com

HC-641-0818