

DESIGNED TO SEAL

# Steri-Tamp®

IV Bag, Vial, & Syringe Seals

- 3X The Adhesive Strength of the Leading Competition
- 100% Sterile Barrier
- Unmistakable Tamper-Evidence



A Brady Business  BRADY

Where positive identification  
meets safe, efficient patient care.

# DESIGNED to Protect Medication from Preparation to Administration

✓ Saves money by allowing return of unused medication to inventory

✓ **Does Not Fall Off**  
Adhesive is 3x stronger than the leading brand<sup>1</sup>, based upon lab testing. Strong adhesive stays sealed, even in cold storage conditions (down to -20° C)

**Stronger Adhesive 3x**

✓ **Tamper-Evident**  
The **ONLY** Tamper-evident seal that features a unique "opened" indicator when removed

✓ **100% Sterile Barrier<sup>2</sup>**

✓ **Single-Use**  
Patented dual-layer design cannot be reapplied

✓ Convenient plastic dispenser box can be wiped down and is USP 797 compliant


✓ Premium medical foil – absolutely no paper in its construction

## It's just a simple seal, right?

If a seal doesn't stay sealed, it has failed its primary purpose. You can't rely upon a seal that works *most* of the time.

That's why Steri-Tamp<sup>®</sup> was invented by a pharmacist who recognized the need for a better seal. Now hospital pharmacies are switching to Steri-Tamp<sup>®</sup> for its greater adhesive strength, reliability, and tamper-evidence.


**Truly Single-Use!**  
**See For Yourself**  
— Pull Tab →



## Bag Port Seals

- Maintains a 100% sterile barrier
- Prevents contamination of medications
- Clearly identifies that IV bag contains added medications to help prevent accidental double dosing
- Provides barrier of protection for pharmacists, technicians, and nurses against exposure to medicine

Product Number	Description	Color	Dimensions	UOM
IVSEALST1	Steri-Tamp® IV Bag Port Seal	Red	1.25" X 2.95"	1000/Roll
IVSEALST2	Steri-Tamp® IV Bag Port Seal	Blue	1.25" X 2.95"	1000/Roll
IVSEALST3	Steri-Tamp® IV Bag Port Seal with "CHEMO" printing	Yellow	1.25" X 2.95"	1000/Roll
IVSEALST4	Steri-Tamp® Belly Button IV Seal for Pfizer/Hospira IV bags	Green	.47" X .78"	1000/Roll



## Vial Seals

- Maintains a 100% sterile barrier
- Prevents contamination of medications
- Convenient, easy-to-use seal for vials
- Tamper-evident feature is highly visible with "OPENED" indicators

Product Number	Description	Color	Dimensions	UOM
VIALSEALST1	Steri-Tamp® Vial Seal, 13MM	Red	.47" X .78"	1000/Roll
VIALSEALST2	Steri-Tamp® Vial Seal, 20MM	Silver	.6" x 1.025"	1000/Roll
VIALSEALST3	Steri-Tamp® Vial Seal, 28MM	Blue	1" x 1.45"	1000/Roll

13mm



20mm



28mm



Get Free Samples & a Quote! Contact Us Today at 800.435.4242  
or visit [www.pdhealthcare.com/steri-tamp](http://www.pdhealthcare.com/steri-tamp)

## Syringe Seals

- Colorful, tamper-evident seal communicates syringe has been sealed by the pharmacy department
- Can also be used as a tamper-evident seal for other containers

**ALLIED**

Product Number	Description	Color	Dimensions	UOM
SYRSEALST1	Steri-Tamp® Syringe Seal	Red	4.25" X .4"	1000/Roll



See Steri-Tamp® application videos online!



## Check Out Other PDC Healthcare Pharmacy Solutions Designed to Help You Improve Patient Safety



### Identify and Manage Medications Across the Continuum of Care with Pharmacy Labeling Solutions

Compatible with pharmacy information systems and automated dispensing systems to provide seamless integration and interoperability with Cerner, MEDITECH, MediWare, OptiFill®, Pyxis® MedStation™, ScriptPro®, VA Hospital, and others.



### Help Prevent Medication Errors with PDC Healthcare's Anesthesia Labels & Tapes

- Prevent medication errors caused by unlabeled or improperly labeled medications
- Widest selection of stock label and tape sizes and formats. Customize any tape or label to meet your facility's specific needs.
- Helps you comply with Joint Commission National Patient Safety Goal NPSG.03.04.01 and ASTM Standard Specification D4774-11

Get Free Samples & a Quote! Contact Us Today at 800.435.4242 or visit [www.pdchealthcare.com/steri-tamp](http://www.pdchealthcare.com/steri-tamp)

**North America & Latin America**  
27770 N. Entertainment Dr., Ste. 200  
Valencia, CA 91355  
Tel 800.435.4242 or 661.257.0233  
Fax 800.321.4409  
[www.pdchealthcare.com](http://www.pdchealthcare.com)  
E-mail [info@pdchealthcare.com](mailto:info@pdchealthcare.com)

**Europe, Middle East & Africa**  
Rue D L'industrie 17  
1400 Nivelles, Belgium  
Tel +32 67 89 56 56  
Fax +32 67 21 77 70  
[www.pdchealthcare.eu](http://www.pdchealthcare.eu)  
E-mail [info@pdchealthcare.eu](mailto:info@pdchealthcare.eu)

**Asia Pacific**  
Tel +1.818.897.1111  
Fax +1 818.686.9317  
E-mail [intl@pdcorp.com](mailto:intl@pdcorp.com)

Where positive identification  
meets safe, efficient patient care.



A Brady Business **BRADY**

HC-476-0515 2,000