

PATIENT EDUCATION

EDUCATIONAL PROCEDURE BROCHURES/BROCHURE HOLDERS

Patient information brochures 2-in-1 English and Spanish in ONE BROCHURE for BIG SAVINGS!

- Improves staff efficiency and instill patient confidence
- Plus, you save money with reduced inventory!
- Imprint area: 2-3/4"W x 1-1/4"H up to 5 lines
- Brochure folds to 3-1/2"W x 8-1/2"H
- 5 pack minimum order
- Typesetting fee will be charged on first or changed order



"Peace of Mind" brochures — What's inside?

"What is CT?"
Describes the procedure

"What are the benefits of CT SCAN?"
Reinforces why the patient needs the procedure

"What are the risks of CT SCAN?"
Educates without alarming the patient

"What will the exam be like?"
Explains what will happen in simple clear language so the patient is less apprehensive

"Illustrations"
Increase understanding with simple drawings

"How should I prepare for the exam?"
Reminds the patient of important preparation steps

"What does the equipment look like?"
Familiarizes the patient to increase their comfort level

"How will I learn the results?"
Confirms the healthcare provider's role in informing the patient

"Remember"
Reinforces key points before and after the procedure

"Questions to ask"
Decreases the likelihood for callbacks with questions


Radiology Brochures			
Product Number	Brochure	Border Color	Qty.
POM-BDE (Non-Imprinted)	Bone Mineral Density Test		100/Pack
POMI-BDE (Imprinted)			
POM-CT (Non-Imprinted)	CT Scan		100/Pack
POMI-CT (Imprinted)			
POM-RAD (Non-Imprinted)	Diagnostic Radiology		100/Pack
POMI-RAD (Imprinted)			
POM-IVP (Non-Imprinted)	Intravenous Pyelogram		100/Pack
POMI-IVP (Imprinted)			
POM-MRI (Non-Imprinted)	MRI		100/Pack
POMI-MRI (Imprinted)			
POM-NUM (Non-Imprinted)	Nuclear Medicine		100/Pack
POMI-NUM (Imprinted)			
POM-PET (Non-Imprinted)	PET - Positive Emission Tomography		100/Pack
POMI-PET (Imprinted)			
POM-ULT (Non-Imprinted)	Ultrasound		100/Pack
POMI-ULT (Imprinted)			
POM-UGI (Non-Imprinted)	Upper GI		100/Pack
POMI-UGI (Imprinted)			

Women's Brochures			
Product Number	Brochure	Border Color	Qty.
POM-BNL (Non-Imprinted)	Breast Needle Localization		100/Pack
POMI-BNL (Imprinted)			
POM-BSE (Non-Imprinted)	Breast Self-Examination		100/Pack
POMI-BSE (Imprinted)			
POM-MAM (Non-Imprinted)	Mammography		100/Pack
POMI-MAM (Imprinted)			
POM-OST (Non-Imprinted)	Osteoporosis		100/Pack
POMI-OST (Imprinted)			
POM-SNB (Non-Imprinted)	Stereotactic Breast Biopsy		100/Pack
POMI-SNB (Imprinted)			
POM-UBB (Non-Imprinted)	Ultrasound-Guided Breast Biopsy		100/Pack
POMI-UBB (Imprinted)			

Call (800) 435-4242 for custom printing

Brochure holders make pertinent information easily accessible to patients!

- Brochure holders are made from clear acrylic
- Available in 8 different sizes with your choice of counter top or wall mount



Product Number	Description	Qty.
PM-C1	1 pocket (Counter Top)	1 Each
PM-W1	1 pocket (Wall Mount)	
PM-C2	2 pocket (Counter Top)	
PM-3-SJ	3 pocket (Counter Top)	
PM-4-SJ	4 pocket (Counter Top)	
PM-C4*	4 pocket (Wall Mount)	
PM-6-SJ	6 pocket (Counter Top)	
PM-W6	6 pocket (Wall Mount)	

Brochures not included with holders. *Clear front, white sides.

Sterile Mini C-Arm Drapes – Sterile

Mini C-Arm Drape fits OEC 6600/6800

- Single piece drape to cover I.I., X-ray Tube and C-Arm

Mini C-Arm Drape (Fluoroscanner/Hologic)

- Single piece drape to cover I.I., X-ray Tube and C-Arm


Product Number	Description	Qty.
CAD2	Mini C-Arm Drape fits OEC 6600/6800 48" x 84"	25/Case
CAD1	Mini C-Arm Drape (Fluoroscanner/Hologic) 52" x 60"	20/Case



Universal C-Arm Drape – Sterile

- Ideal for cases where only the region above the table needs to be draped
- Designed to fit nearly all C-Arm designs
- Clear with poly straps
- One piece design

Product Number	Description	Qty.
CAD3	42" x 74"	20/Case



X-ray Tube Covers and Universal Halfbag Covers – Sterile

9" and 12" X-ray Tube Covers

- 9" Band Bag (30" x 32")
- 12" Band Bag (36" x 36")

9" Universal Halfbag Cover

- Fits the I.I. and midway down the C-Arm

Product Number	Description	Qty.
CABG1	9" Band Bag (30" x 32")	25/Case
CABG2	12" Band Bag (36" x 36")	25/Case
CABG3	9" Halfbag Cover (36" x 80")	20/Case




C-Arm Kit – Sterile

Kit Includes:

- C-Arm drape with 5 mounting clips
- 36" x 36" I.I. cover
- 36" x 36" X-ray tube cover
- 24" x 20" Footswitch cover

Product Number	Description	Qty.
CA12	C-Arm Kit	20 Kits/Case



Microscope Drapes – Sterile

- Made with strong, drapable polyethylene film with monocular or stereo observer tubes
- Include straps to gather excess poly, and re-attachable fiber-free ocular closures

Product Number	Accommodates Microscope	Fits Lens Size	Size	Qty.
MSD12	Zeiss OPMI	65mm	42" x 120"	10/Case
MSD42	Zeiss MD	65mm	42" x 120"	
MSD44	Zeiss MD	65mm	44" x 120"	
MSD54	Leica Wild	68mm	54" x 150"	
MSD65	Zeiss MD	65mm	44" x 65"	



PATIENT/CAREGIVER SAFETY

PROBE COVERS

General Purpose Probe Covers

Reduce cross-contamination risk when performing routine puncture and drainage procedures. Select from flat fold or telescopic fold, CIV-Flex or Polyethylene.



Intraoperative Probe Covers

Extend the sterile field in operating rooms and surgi-suites by covering surgical transducers. Accordion fold.



Flat Fold Endocavity Probe Covers

Use in all prostate, rectal and vaginal procedures. Provides a barrier to reduce cross-contamination risks. Select from flat fold or rolled.



Rolled Endocavity Probe Covers

These CIV-Flex or Polyethylene probe covers are acoustically transparent, soft and comfortable for patients.



Product Number	Cross Reference	Function	Packaging	Contents	Material	Fold	End	Dimensions (CM)	Qty.
PRB-517	610-002	General Purpose	S, IP	gel-pkt, w/bands	CIV-Flex	Flat Fold	3D end	14 x 61	24/Box
PRB-542	610-542	General Purpose	S, IP	gel-pkt, w/bands	CIV-Flex	Telescopic Fold	3D end	14 x 91.5	24/Box
PRB-836	610-836	Endocavity	NS, IP	w/bands	Neo-Guard	Rolled	—	2.6 x 30	24/Box
PRB-843	610-843	Endocavity	S, IP	w/bands	Neo-Guard	Rolled	—	2.6 x 30	24/Box
PRB-838	610-838	Endocavity	NS, IP	w/bands	Neo-Guard	Rolled	—	4 x 30	24/Box
PRB-844	610-844	Endocavity	S, IP	w/bands	Neo-Guard	Rolled	—	4 x 30	24/Box
PRB-521	610-007	Endocavity	NS, IP	w/bands	CIV-Flex	Flat Fold	—	11.9 > 4.6 taper x 61	24/Box
PRB-610	610-792	Endocavity	NS, IP	Pre-Gelled	Polyethylene	Flat Fold	—	6.3 > 4.6 taper x 24.1	100/Box
PRB-367	610-367	Intraoperative	S, IP	gel-pkt, w/bands	Polyethylene	Accordion Fold	—	15.2 > 7.6 taper x 244	12/Box
PRB-689	610-689	BANDS	NS, BULK	BANDS (orange)	Kraton	—	—	5/8"ID (inner diameter)	100/Box

S = Sterile Pack, NS = Non-Sterile Pack, IP = Individually Packaged

FLAT FOLD - Folded end over end

TELESCOPIC FOLD - Material is folded within itself, but not rolled. This facilitates both gel application, and transducer insertion.

ACCORDION FOLD - Gathered, the Surgi Intraoperative covers provide distortion-free scanning and extend the length of the transducer cord for up to eight feet of sterile protection in the OR suite.



For our complete line of Equipment Covers, visit us at www.pdhealthcare.com

DID YOU KNOW?

According to the CDC, 1.7 million people a year contract hospital-related infection, resulting in 99,000 deaths.

Source: www.cdc.gov/HAI/pdfs/hai/infections_deaths.pdf

SANI-CLOTH® PLUS

low alcohol formula

- Kills bacteria and viruses
- Effective against 15 micro-organisms including TB, HBV and MRSA
- Contact time is 5 minutes or less

Product Number	Description	Wipe Size	Qty.
INF-635	Sani-Cloth® Plus	6-1/2" x 6-1/4"	160 Wipes/Tub



SUPER SANI-CLOTH®

high alcohol formula

- Kills bacteria, viruses and yeast
- Effective against 26 micro-organisms including TB, HBV and MRSA
- Contact time is 5 minutes or less

Product Number	Description	Wipe Size	Qty.
INF-639	Super Sani-Cloth®	6" x 6-3/4"	160 Wipes/Tub
INF-636	Super Sani-Cloth®	11-1/2" x 11-3/4"	50 Packets/Box



SANI-CLOTH® HB

alcohol-free formula

- Kills bacteria, viruses and fungi
- Effective against over 100 micro-organisms including TB, HBV and MRSA
- Contact time is 10 minutes or less

Product Number	Description	Wipe Size	Qty.
INFH635	Sani-Cloth® HB	6" x 6-3/4"	160 Wipes/Tub
INFH637*	Sani-Cloth® HB	8" x 14"	65 Wipes/Tub

*Not shown



CaviCide® Surface Cleaner and Wipes

- Kills TB in 3 minutes and MRSA, HIV-1, HBV and HCV in 2 minutes
- Kills C-DIFF vegetative cells in 2 minutes
- EPA approved to be effective against H1N1
- Low alcohol (17.20% Isopropanol)
- Available in wipes or spray

Product Number	Description	Wipe Size	Qty.
INF311	CaviWipes®	6" x 6-3/4"	160 Wipes/Tub
INF312*	CaviWipes®	9" x 12"	65 Wipes/Tub
UPI-43	CaviCide® Spray	24 oz.	1 Bottle/Each

*Not shown



T-Spray II™ Transducer Detergent Disinfectant Spray

- Requires no scrubbing – simply wipe clean
- Leaves no grit or soap scum
- Will not cause swelling of transducer membrane
- Contains no phenol, alcohol or glutaraldehyde
- Great for mammography

Product Number	Size	Qty
UPI-39	250 ml. Bottle	1 Each
UPI-37*	250 ml. Bottle	12/Case

*Not shown



PATIENT/CAREGIVER SAFETY

BIOHAZARD BAGS/TABLE PADS/PATIENT TRANSFER BOARD/SANDBAGS



More sizes, closures and styles available!

Ziplock with paperwork pouch


BIOBG8

Specimen Transport Bags


- Pre-printed with "BIOHAZARD" for easy identification
- Paperwork pouch designed with a flap cover to secure paperwork during transfer
- 2 mil thickness
- More sizes, closures and styles available

Product Number	Description	Size	Qty.
BIOBG8	Ziplock Closure	8" x 8"	1,000/Case
BIOBG9	Ziplock Closure	6" x 9"	

Call (800) 435-4242 for custom printing



2152




Use this guide to specify your dimensions.


Standard Foam Table Pads

- Durable almond colored vinyl covering for strength and comfort
- Convenient zipper is included for easy removal and cleaning
- Optional hook and loop or velcro closure available

Product Number	Description	Size	Qty.
2152	Standard Table Pads	1"H x 24"W x 72"L	1 Each
2153	Standard Table Pads	1"H x 27"W x 78"L	
2155	Standard Table Pads	1-1/2"H x 27"W x 78"L	



STJ-940



STJ940R

Handy storage rack for maximum convenience!

Patient Transfer Board & Storage Rack

- Anti-static, stress relieved polyethylene plastic for easy movement
- Radiolucent so that you can X-ray your patient while still on the board
- Wall mount storage rack available

Product Number	Description	Size	Qty.
STJ-940	Patient Transfer Board	23"W x 70"L x 1/4"H	1 Each
STJ940R	Wood Storage Rack for STJ-940	—	



Heavy-duty nylon straps are softer and gentler on patient's hands

SB-750

Sandbag Set*:

- Two, 10 lb. 11" x 11" Sandbags with Handles
- Two, 5 lb. 5" x 16" Sandbags with Eyelet
- Two, 2-1/2 lb. 6-1/2" x 9" Sandbags with Eyelet

Available Colors:



"Soft Strap" sand bags keep patients stationary.

- Heat sealed inner sleeve to prevent sand leakage
- Heavy gauge reinforced vinyl outer cover
- Sewn seams are tear resistant
- Straps are long enough to go around the wrist

Product Number	Description	Color	Size	Qty
SB-750	Sandbag Set (6 pc. set)	Black	*Assorted	1 Set
SB-736	5 lb. Sandbags with Handles (2 per set)	Black	9" x 9"	1 Pair
SB-738	5 lb. Sandbags with Handles (2 per set)	Green	9" x 9"	1 Pair
SB-739	5 lb. Sandbags with Handles (2 per set)	Blue	9" x 9"	1 Pair
SB-726	10 lb. Sandbags with Handles (2 per set)	Black	11" x 11"	1 Pair
SB-727	10 lb. Sandbags with Handles (2 per set)	Red	11" x 11"	1 Pair
SB-729	10 lb. Sandbags with Handles (2 per set)	Blue	11" x 11"	1 Pair