

# AUTOCLAVE LABELS

## Increase Workflow Efficiency

Autoclave labels enable the proper identification of equipment, instruments and supplies as they go through the autoclave process. They are designed to be applied onto a tray or tag to identify who owns the surgical equipment and instruments and their intended procedures. This helps to improve communication and efficiency, getting the items into surgery more quickly.

- Durable and lab-tested to withstand autoclave and common disinfectant solutions
- Tested to withstand vaporized hydration peroxide (VHP) sterilization, including STERRAD® Sterilization Systems\*
- Made of premium synthetic material for thermal transfer printing
- Withstand temperatures up to 300°F
- Remain adhered through autoclave without peeling or falling off

\*STERRAD® is a registered trademark owned by Advanced Sterilization Products (ASP)

### Thermal Transfer Autoclave Labels

SKU NUMBER	LABEL DIMENSIONS	CORE SIZE
TTATCLV2	2" x 1"	1"
TTATCLV23	2" x 1"	3"
TTATCLV21	2" x 1-¼"	1"
TTATCLV213	2" x 1-¼"	3"
TTATCLV25	2" x 1-½"	1"
TTATCLV31	3" x 1"	1"
TTATCLV313	3" x 1"	3"
TTATCLV32	3" x 2"	1"
TTATCLV323	3" x 2"	3"
TTATCLV4112	4" x 1-½"	1"
TTATCLV41123	4" x 1-½"	3"
TTATCLV42	4" x 2"	1"

### Thermal Transfer Printer Ribbons


SKU NUMBER	RIBBON SIZE	CORE SIZE	UNIT OF MEASURE
TTR-FED110G30	4.33" x 243'	1/2"	36 ROLLS/BOX
TTR-FRD11030-12	4.33" x 243'	1"	12 ROLLS/BOX
TTR-FED0644S	2.5" x 243'	1/2"	12 ROLLS/BOX
TTR-FED110G30	4.33" x 243'	1/2"	36 ROLLS/BOX

Need help identifying the label and ribbon compatible with your printer? Talk to a PDC representative today.



Order Kit # TTATCLV31-KIT to receive receive both labels (TTATCLV31) and the compatible ribbons (TTR-FED110G30)

# OTHER STERILE PROCESSING SOLUTIONS



Duraholder® Surgical Instrument Holders



Ident-Alert® Tray & Scope Tags



Autoclave Tape



Communication Labels



Handheld Label Dispensers & Labels

Learn more  
**800.435.4242 | [pdhealthcare.com](http://pdhealthcare.com)**

